

Product Data Sheet

DIAION™ HP20

DIAION™ HP20 is based on a unique rigid polystyrene/divinylbenzene matrix. A controlled pore size distribution and large surface area offer excellent resolution and the capacity for a wide range of molecules, from small peptides and oligonucleotides up to large proteins.

DIAION™ HP20 is characterized by:

- >> Unique pore size distribution
- >> Excellent batch-to-batch reproducibility
- >> Wide application
- >> High chemical and physical stability
- >> Excellent pressure/flow characteristics


Physical and chemical properties

Grade Name	DIAION™ HP20	
Bead Form	Spherical, porous	
Matrix	Polystyrene/divinylbenzene	
Chemical Structure	$\begin{array}{c} \text{---CH}_2\text{---CH---CH}_2\text{---CH---} \\ \qquad \qquad \\ \text{C}_6\text{H}_5 \qquad \text{C}_6\text{H}_5 \\ \qquad \qquad \\ \text{---CH---CH}_2\text{---} \end{array}$	
Shipping Density*	g/L	690
Water Content	%	55 - 65
Particle Size Distribution thr. 250 μm	%	10 max.
Effective Size	mm	0.25 min.
Uniformity Coefficient	-	1.6 max.
Particle Density*	g/mL	1.01
Specific Surface Area*	m ² /g	590
Pore Volume*	mL/g	1.3
Pore Radius*	Å	290

Note : properties with a mark "*" are referential data.

Swelling ratio in various solvents

Methanol	1.13
Ethanol	1.24
2-Propanol	1.17
Acetone	1.24
Toluene	1.25
Acetonitrile	1.17
Water	1.00

Pore size distribution**Fig. 1 Pore size distribution of HP20****Recommended Operating Conditions**

Maximum Operating Temperature	°C	130
Operating pH Range		0 - 14
Minimum Bed Depth	mm	800
Flow rate	BV/h	Loading 0.5 - 5
	BV/h	Displacement 0.5 - 2
	BV/h	Regeneration 0.5 - 2
	BV/h	Rince 1 - 5

Regenerant

- Organic solvents for hydrophobic compounds
- Bases for acidic compounds
- Acids for basic compounds
- Buffer solution for pH sensitive compounds
- Water for an ionic solution
- Hot steam for volatile compounds

Hydraulic Characteristics

The approximate pressure drop at various temperatures and flow rates for each meter of bed depth of DIAION™ HP20 resin in normal down flow operation is shown in the graphs below.




Fig. 2 Pressure Drop of HP20

FDA status

DIAION™ HP20 may be used to process food and beverage products and isolate specialized food additives as intended and such used may be said to fully comply with the Federal Food, Drug, and Cosmetic Act.

Applications

- Purification of small peptides, oligonucleotides and proteins
- Adsorption of vitamins, antibiotics, enzymes, steroids and other substance from fermentation solutions
- Decolorization of various sugar solutions
- Adsorption of fatty acids
- Removal of phenol
- Adsorption of various perfume
- Decolorization and purification of various chemicals

Storage condition

Synthetic adsorbents are recommended to store properly in order to avoid a high risk for mold growth. The proper storage typically allows any synthetic adsorbent resin to last for a year after production before onset of any such growth.

The best storage condition is with 20% of alcohol such as ethanol or isopropanol. A 10% or higher concentration of salt solution, such as NaCl, is also recommended to preserve new or used resin for long storage.

In case salt cannot be used, a 0.01 to 0.02 N of NaOH solution could be accepted as mold cannot withstand survival at pH higher than 12.

Storage at freezing temperature should be avoided at all cost as it may cause breakage or crush of resin particles.

Notice

This information are given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. The application, use and processing of our products are beyond our control and therefore your own responsibility.