

Product Data Sheet


SEPABEADS™ SP850

SEPABEADS™ SP850 is highly porous styrenic adsorbents. It has much larger surface area and a narrower pore size distribution than DIAION™ HP20. It has also smaller pore radius than SEPABEADS™ SP825L. It offers higher capacity for small molecules. This grade is recommended for adsorption, desalting and decolorization.

SEPABEADS™ SP850 is characterized by:

- >> Unique pore size distribution
- >> Excellent batch-to-batch reproducibly
- >> Wide application
- >> High chemical and physical stability
- >> Excellent pressure/flow characteristics

Physical and chemical properties

| | | |
|--|--|-----------|
| Grade Name | SEPABEADS™ SP850 | |
| Bead form | Spherical, porous | |
| Matrix | Polystyrene/divinylbenzene | |
| Chemical Structure |  | |
| Whole beads count | - | 95 min. |
| Shipping Density* | g/L | 690 |
| Water content | % | 46 - 52 |
| Particle Size Distribution thr. 250 µm | % | 10 max. |
| Effective size | mm | 0.25 min. |
| Uniformity Coefficient | - | 1.6 max. |
| Particle Density* | g/mL | 1.01 |
| Specific Surface Area* | m ² /g | 930 |
| Pore Volume* | mL/g | 1.1 |
| Pore Radius* | Å | 45 |

Note : properties with a mark "*" are referential data.

Swelling ratio in various solvents

| | |
|--------------|------|
| Methanol | 1.15 |
| Ethanol | 1.17 |
| 2-Propanol | 1.19 |
| Acetone | 1.17 |
| Toluene | 1.15 |
| Acetonitrile | 1.15 |
| Water | 1.00 |

Pore size distribution




Fig. 1 Pore size distribution of SP850

Recommended Operating Conditions

| | | |
|-------------------------------|------|----------------------|
| Maximum Operating Temperature | °C | 130 |
| Operating pH Range | | 0 - 14 |
| Minimum Bed Depth | mm | 800 |
| Flow rate | BV/h | Loading 0.5 - 5 |
| | BV/h | Displacement 0.5 - 2 |
| | BV/h | Regeneration 0.5 - 2 |
| | BV/h | Rince 1 - 5 |

Regenerant

Organic solvents for hydrophobic compounds

Bases for acidic compounds

Acids for basic compounds

Buffer solution for pH sensitive compounds

Water for an ionic solution

Hot steam for volatile compounds

Hydraulic Characteristics

The approximate pressure drop at various temperatures and flow rates for each meter of bed depth of SEPABEADS™ SP850 resin in normal down flow operation is shown in the graph below.




Fig. 2 Pressure Drop of SP850

Applications

- Purification of Cephalosporin C
- Purification of small peptides, oligonucleotides and proteins
- Adsorption of vitamins, antibiotics, enzymes, steroids and other substance from fermentation solutions
- Decolorization of various sugar solutions
- Adsorption of fatty acids
- Removal of phenol
- Adsorption of various perfume
- Decolorization and purification of various chemicals

Storage condition

Synthetic adsorbents are recommended to store properly in order to avoid a high risk for mold growth. The proper storage typically allows any synthetic adsorbent resin to last for a year after production before onset of any such growth.

The best storage condition is with 20% of alcohol such as ethanol or isopropanol. A 10% or higher concentration of salt solution, such as NaCl, is also recommended to preserve new or used resin for long storage.

In case salt cannot be used, a 0.01 to 0.02 N of NaOH solution could be accepted as mold cannot withstand survival at pH higher than 12.

Storage at freezing temperature should be avoided at all cost as it may cause breakage or crush of resin particles.

Notice

This information are given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. The application, use and processing of our products are beyond our control and therefore your own responsibility.